




Grease-free distributor for vacuum and inert gas



Application: Laboratories with oxygen- and air sensitive substances, with ignitable substances and pollution-sensitive products (e.g. pollution due to grinding grease)

- Advantages:
- Grease free
 - Highest vacuum tightness due to
 - High vacuum flanges
 - NORMAG spindle valves
 - High operating safety due to
 - different coloured valve caps for vacuum and inert gas
 - operator-convenient design
 - High operating safety due to
 - use of sophisticated DURAN-glass
 - NORMAG overflow valve with defined actuation pressure and blow-off capacity
 - Customer-specific design with regard to
 - connections
 - sizing / adaptation to local conditions